



Savanna® - Respiratory Viral Panel 4 (RVP-4)

Project Status and next steps

- Dr. rer. nat. Jürgen Becker
- Scientific Marketing and Clinical Value Manager EMEA

Content: Results, status and next steps related to

- (I) Clinical Performance Studies
- (II) Cost benefit Evaluation Savanna RVP-4
- (IV) Publication plan and Next steps



(I) Completed Clinical Performance Studies - OVERVIEW

Site	Country	Investigator	Comparator	Periode
multiple	USA	multiple	BioFire Respiratory Panel 2.1	2021
LKH Feldkrich	Austria	Dr. Andreas Bösl	Roche Cobas SRAS CoV-2	Apr. 2022
Pievesistina	Italy	Pr. Vittorio Sambri	Seegene Allplex	July 2022
			Abbott Alinity SARS CoV-2	
			Hologic Aptima SARS CoV-2 TMA	
CHU Toulouse	France	Pr. Jaques Izopet	Hologic Influenza A, B RSV on Panther	Dec. 2022 and Jan. 2023
			Hologic SARS CoV-2, Flu A, Flu B, RSV on Panther	
			Hologic Aptima SARS CoV-2 TMA	2025
Chur	Swisse	Dr. Riccarda Capaul	Xpert Xpress SARS CoV-2, Flu A, Flu B, RSV	Jan and Feb. 2023

(I) Completed Clinical Performance Studies - RESULTS

Ctudy	Positive % Agreement		negative % Agreement		Overall %-	N		
Study	Value	95% CI	Value	95% CI	Agreement	all	Pos	Neg
USA (5 Centers)	99,44%	96,43 – 99,97	100%	99,34 – 100	99,89%	896	178	718
Feldkirch (Austria)	93.62%	87.87 – 96.85	97.98%	94.57 – 99.35	96,17%	339	141	198
Pievesistina (Italy)	93,83%	85.55 – 97.71	98,29%	93.35 – 99.70	96,46%	198	81	117
Toulouse (France)	92,11%	87,29 – 95,28	99,81%	99,23 – 99,97	98,56%	1251	203	1048
Chur (Swiss)	92,59%	81.26 – 97.60	100%	98.98 – 100	99,23%	520	54	466
All (1-5)	94,67%	92.59 – 96.21	99,69%	99.36 – 99,85	98,66%	3204	657	2547

(II) Cost Benetit-Evaluation Savanna RVP-4







Article

Cost–Benefit of Real-Time Multiplex PCR Testing of SARS-CoV-2 in German Hospitals

Roland Diel 1,2,3,* and Albert Nienhaus 3,4





Manuscript statistic
- 540 views in 4 weeks

Results:

In probabilistic sensitivity analysis, implementing the Savanna® test saved, on average, €107 as compared to applying the clinical-judgement-only strategy. A revenue loss of €735 can be avoided when SARS-CoV-2 infection in patients coming unplanned to the hospital due to other acute illnesses are excluded immediately by POCT.

Conclusions:

Using highly sensitive and specific PCR-POCT in patients suspected of COVID-19 infection at German ERs may significantly reduce hospital expenditures.





ENDE